



DeGen
Medical

Cyclops™ **Anterior Cervical** **Plate System**

The Cyclops™ **Anterior Cervical Plate System**

provides an exceptionally short, low profile plate with hyper-angulated screws for use in anterior cervical discectomy and fusion.

- Large window enables visualization of the end plate for easy plate selection
- 1-5 level plates in multiple lengths
- Single level plates available in one-millimeter increments from 8mm - 13mm
- 2.2mm thin plate reduces the risk of dysphagia
- Convenient central cam design provides a lock for both screws at each level

DeGenMedical.com



DeGen
Medical

