



## <sup>3</sup>e MIS™ Pedicle Screw System



### The <sup>3</sup>e MIS™ Pedicle Screw System

provides multiple head and screw options that allow for a customized approach to each case. The modular heads and screws of this system are compatible with modular heads and screws of the F1 MPS™ Modular Pedicle Screw System.

- 56° multi-axial angulation
- Double lead thread for cortical and cancellous bone purchase
- Cannulated diameters from 4.5mm to 10.0mm in 0.5mm increments
- Solid diameters from 4.0mm to 10.0mm in 0.5mm increments
- Self-tapping flute

[DeGenMedical.com](http://DeGenMedical.com)



**DeGen**  
Medical