




DeGen
Medical

Latitude-C™ Cervical Interbody Spacer System

The Latitude-C™ Cervical Interbody Spacer System

is designed as a laterally-angled wall interbody fusion device. It has the option of five different footprints with different angles for uncinete processes and up to ten different heights which allows the surgeon to choose among a variety of anatomic shapes for every cervical level and patient.

- Restoration of disc height and lordosis
- Compliments anatomical shape
- Large graft volume
- Size options:
 - 12 x 12
 - 14 x 12
 - 16 x 14
 - 18 x 14
 - 20 x 15



DeGenMedical.com



DeGen
Medical

