



DeGen
Medical

Impulse AMTM

Lumbar Interbody Fusion System

THE IMPULSE AMTM INTERBODY FUSION SYSTEM

offers an interbody spacer that is additively manufactured from DeGen Medical's proprietary titanium, Puri-TiTM.


These Implants are offered in a wide range of geometries, footprints, and heights to match complex patient anatomy. The Diamond lattice on the lateral walls allow for increased visualization while maintaining mechanical integrity.

THE IMPULSE AMTM INTERBODY SPACER


| WIDTHS ¹ : | LENGTHS: | HEIGHTS: | GEOMETRIES: |
|-----------------------|----------|----------------------------------|--------------|
| 8mm | 22mm | 7mm - 16mm | Neutral |
| 10mm | 26mm | (in 1mm increments) ² | Anatomic |
| 12mm | 32mm | | Lordotic 8° |
| | 36mm | | Lordotic 12° |

¹Widths are identified as whole numbers

²32 and 36 long, 12 degrees excluded 7,8 mm heights



Top Level



Bottom Level

